


LabChart
LIGHTNING

SDK Documentation

LabChart
LIGHTNING

Overview (single device)



The **DeviceClass** is implemented by your plugin code. Each device class object can support a single device, or a range of related devices that can handle the same types of settings.

On launch, LabChart Lightning loads all plugins.

When LabChart Lightning starts, it calls **getDeviceClasses()** on your plugin. Each plugin can, if desired, implement more than one **DeviceClass**. Note that Lightning will not let settings from one **DeviceClass** be applied to a device from a different class (this is essentially the defining feature of a **DeviceClass**).

Typically there would only be one **DeviceClass** per plugin (unless the plugin needs to support devices of very different types).

This must return an object array which Lightning can query by calling **DeviceClass.checkDeviceIsPresent()** for each possible device connection.


A **PhysicalDevice** object is created by your plugin code and represents an instance of a single, connected hardware device. Your plugin creates **PhysicalDevice** objects for as many devices that are discovered.

Multiple LabChart Lightning recordings can use the same **PhysicalDevice**, but only one can sample with that device at any time.

Each recording in Lightning holds the settings for that device for that recording and allows the getting and setting of device settings specific to the recording.

Lightning uses the **DeviceClass** supplied by the plugin to create these **DeviceProxy** objects as required.

Overview (multiple devices)



The **DeviceClass** is implemented by your plugin code. Each device class object can support a single device, or a range of related devices that can handle the same types of settings.

On launch, LabChart Lightning loads all plugins.

When LabChart Lightning starts, it calls **getDeviceClasses()** on your plugin. Each plugin can, if desired, implement more than one **DeviceClass**. Note that Lightning will not let settings from one **DeviceClass** be applied to a device from a different class (this is essentially the defining feature of a **DeviceClass**).

Typically there would only be one **DeviceClass** per plugin (unless the plugin needs to support devices of very different types).

This must return an object array which Lightning can query by calling **DeviceClass.checkDeviceIsPresent()** for each possible device connection.

A **PhysicalDevice** object is created by your plugin code and represents an instance of a single, connected hardware device. Your plugin creates **PhysicalDevice** objects for as many devices that are discovered.

Multiple LabChart Lightning recordings can use the same **PhysicalDevice**, but only one can sample with that device at any time.

Each recording in Lightning holds the settings for that device for that recording and allows the getting and setting of device settings specific to the recording.

Lightning uses the **DeviceClass** supplied by the plugin to create these **DeviceProxy** objects as required.

Plugin location and structure

Location

C:\Users\[USERNAME]\Documents\LabChart Lightning\Plugins\Devices

Windows

~/Documents/LabChart Lightning/Plugins/Devices

macOS

Structure/contents


 Plugins

 Devices

 plugin-name

 plugin-name.js **(required)**

The entry point for your plugin's code. Split up your code into multiple files and import/reference them in this file.

 Additional files **(optional)**

Split up your code into as many (or as few) files as necessary. Import/reference these files in the *plugin-name.js* file.

What your plugin needs to implement

You will need to create the following TypeScript classes as well as their required methods.

DeviceClass

(extends `IDeviceClass`)

The **DeviceClass** object represents this set of devices and can find and create **PhysicalDevice** objects of its class, as well as the **ProxyDevice** objects.

All devices in this class can understand the same settings. An important function of the **DeviceClass** is enable Lightning to match hardware devices to previous device settings (i.e. **ProxyDevices**) when reopening recordings by matching properties such as serial numbers and device capabilities.

The **DeviceClass** also contains properties shared by all devices in the class.

ProxyDevice

(extends `IProxyDevice`)

The **ProxyDevice** object is created for each recording and manages the device settings and sampling for its recording.

PhysicalDevice

(extends `OpenPhysicalDevice`)

PhysicalDevice is a representation of the connected hardware device. Once created, **PhysicalDevices** last for the entire Lightning session.

The **PhysicalDevice** holds onto the **Parser** and holds the capabilities of its hardware device, e.g. maximum number of inputs, the serial number and other properties that can be used to identify that particular piece of hardware for matching settings to devices when recordings are reopened.

What your plugin needs to implement

Inputs vs streams

Input

Each Input corresponds to a hardware input on the device (e.g. a BNC connector on a PowerLab), so it has a Range, typically in V, mV or uV which depends on the gain of the input (which may be fixed or adjustable).

An input may have other device dependent properties such as AC coupling, Differential (vs single ended), etc, which are stored in the **InputSettings** object.

Stream

Each **ProxyDevice** provides sample data to its recording via one or more data streams. The number of streams can depend both on the capabilities of the hardware and the settings.

The number of streams need not equal the number of hardware inputs, e.g. the device may calculate more than computed data stream from each input.

Required classes:

DeviceClass (IDeviceClass)

The **DeviceClass** object represents a set of devices that can share the same types of settings and can find and create **PhysicalDevice** objects of its class, as well as the **ProxyDevice** objects.

ProxyDevice (IProxyDevice)

The **ProxyDevice** object is created for each recording that uses the device and manages hardware settings and sampling and access to the **PhysicalDevice** (only one recording can access the **PhysicalDevice** at any one time).

PhysicalDevice (OpenPhysicalDevice)

PhysicalDevice is a representation of the connected hardware device and owns the parser and the device connection.

Parser

A custom object, designed by you, that handles parsing of data returned from the device and adds the resulting samples to the device's output streams (signals).

Other classes:

InputSettings

The input corresponds to a hardware input on the device (e.g. a BNC connector on a PowerLab), so it has a Range, typically in V, mV or uV which depends on the gain of the input (which may be fixed or adjustable).

An input may have other device dependent properties such as AC coupling, Differential (vs single ended), etc.

Sometimes a device provide more streams (signals) than it has inputs, generally because it calculates multiple signals from one or more of its inputs.

Often these computed streams have different units or sample rates from the input(s).

StreamSettings (IDeviceStreamApi)

The settings for each Stream for each recording are stored in a **StreamSettings** object. The **ProxyDevice** holds an array of these.

A **StreamSettings** object can hold arbitrary, device or stream specific settings in addition to the standard ones: "enabled", "samplesPerSec" and "inputSettings" (which is of type **InputSettings**).

DeviceStreamConfigurationImpl (IDeviceStreamConfiguration)

This object holds the information Quark needs to know about the units and data format of the Stream. The **unitsInfo** will change if the gain of the stream's input changes.



DeviceClass

(extends `IDeviceClass`)

The `DeviceClass` object represents this set of devices and can find and create `PhysicalDevice` objects of its class, as well as the `ProxyDevice` objects.

Identifying the device

`[]` `getClassId(): string`

Should return a universally unique identifier (uuid) specific to this device class.

`[]` `getDeviceClassName(): string`

Should return a human readable name for the `DeviceClass` to be shown in Lightning UI.

`[]` `getDeviceConnectionType(): TDeviceConnectionType;`

Specifies the type of the connections Quark should create and pass to `checkDeviceIsPresent()` to support device discovery. E.g. `kDevConTypeSerialPort` or `kDevConTypeSerialOverBluetooth`.

Handling connections

`[]` `createProxyDevice(
 quarkProxy: ProxyDeviceSys | null,
 physicalDevice: OpenPhysicalDevice | null
): IProxyDevice`

`[]` `checkDeviceIsPresent(
 deviceConnection: DuplexDeviceConnection,
 callback: (error: Error | null, device: OpenPhysicalDevice | null) => void
): void`

`[]` `indexOfBestMatchingDevice(
 descriptor: OpenPhysicalDeviceDescriptor,
 availablePhysDevices: OpenPhysicalDeviceDescriptor[]
): number`

Called when deciding (typically when loading a recording) which physical device should be used for the specified device proxy in a recording.

This function should run through all the available `PhysicalDevices` comparing them against the capabilities and device id properties passed in the descriptor parameter in order to choose the best match.

If a valid index is returned, that `PhysicalDevice` will be removed from the available `PhysicalDevices` passed in future calls.

Handling errors

`[]` `onError(err: Error): void` optional

Called when a sampling or connection error is detected. You can choose to ignore these.

Other

`[]` `release?(): void` optional

Called when the app shuts down. Chance to release any resources acquired during this object's lifetime.



PhysicalDevice

(extends `OpenPhysicalDevice`)

PhysicalDevice is a representation of the connected hardware device.

PhysicalDevices need to own a custom parser object which typically holds onto the **PhysicalDevice's DeviceConnection** so that the parser can always receive any bytes arriving from the device.

What you must implement

① `getDeviceName(): string;`

① `getNumberOfAnalogInputs(): number;`

① `getNumberOfAnalogStreams(): number;`

Commonly, this will be the same as the number of analog inputs, unless the device can generate multiple calculated data streams from its hardware inputs.

① `getDescriptor(): OpenPhysicalDeviceDescriptor;`

① `release(): void;`

Since the **PhysicalDevice** owns its **DeviceConnection**, `release()` should call `onStreamDestroy()` and then `release()` on the **deviceConnection**, (which is normally held onto by the **PhysicalDevice's** parser).



ProxyDevice

(extends IProxyDevice)

The **ProxyDevice** object is created for each recording and manages access to the **PhysicalDevice** which may be shared by multiple recordings, each having different settings.

What you must implement

`[] getOutBufferInputIndices(): Int32Array;`

`[] getDeviceName(): string;`

`[] getNumberOfAnalogStreams(): number;`

`[] getLastError(): string;`

`[] setPhysicalDevice(physicalDevice: OpenPhysicalDevice): boolean;`

`[] connectToPhysicalDevice(): boolean;`

Called from Quark to allow this proxy to communicate with the device.

`[] disconnectFromPhysicalDevice(): void;`

Called from Quark to stop this proxy communicating with the device to allow another proxy to use the device.

`[] prepareForSampling(bufferSizeInSecs: number): boolean;`

Allocate StreamRingBuffers buffers.

`[] startSampling(): boolean;`

`[] onSamplingStarted(): void;`

`[] onSamplingUpdate(): void;`

`[] stopSampling(): boolean;`

`[] onSamplingStopped(errorMsg: string): void;`

`[] cleanupAfterSampling(): boolean;`

Release buffers and cleanup.

`P isSampling: boolean;`


`P outputStreamBuffers: StreamRingBuffer[];`

Inputs, streams, and signals.

Inputs

Each Input corresponds to a hardware input on the device (e.g. a BNC connector on a PowerLab), so it has a Range, typically in V, mV or μV which depends on the gain of the input (which may be fixed or adjustable).


An input may have other device dependent properties such as AC coupling, Differential (vs single ended), etc, which are stored in the InputSettings object.



Streams


Each ProxyDevice provides sample data to its recording via one or more data streams. The number of streams can depend both on the capabilities of the hardware and the settings.

The number of streams need not equal the number of hardware inputs, e.g. the device may calculate more than computed data stream from each input.



Signals

Lightning maps the data streams to signals and exposes these signals to users.



The connection process



Lightning is launched.

On launch, LabChart Lightning loads all plugins by calling `getDeviceClasses()` on each plugin. This returns an instance of the plugin `DeviceClass` implementation. e.g. `[return new DeviceClass()]`;



Lightning scans for hardware.

For each device connection that is not already being used by a physical device:

`checkDeviceIsPresent(deviceConnection, callback)` is called on the `DeviceClass` implementation.

It is up to the `checkDeviceIsPresent` function within the script to decide whether its hardware is present, e.g. the device might use a generic serial to USB chip such as FTDI such that the vendor, product and manufacturer information (vid, pid etc) are not sufficient to identify the actual device type, in which case it will need to try communicating with it.

If a device of the right type is detected, a new `PhysicalDevice` object is instantiated and returned to LabChart Lightning via callback.



DeviceClass
Acme devices plugin


PhysicalDevice objects created.

When physical devices are created during search for devices, any recordings using devices of that class are notified so that they can connect any proxies of that class to the physical devices.

When recordings are created, one of two things happen:

- (1) If a recording is created with default settings, a `DeviceProxy` is created for the default device for use by the recording.
- (2) If a recording is cloned from an existing recording, the `DeviceProxies` within the existing recording are cloned for use by the new recording.


The Sampling Process



Places your plugin can present UI


Configuring signals

Each signal from your device can be configured using the Signal properties popover in the ChartView. Your plugin script can define custom UI for this popover.




Device connection

Devices handled by your plugin script will be shown in either the Add Device list (if it is not used by the active recording) or the Recording Devices list (if it is being used). Name, status, and button to link to custom user interface (if any exists) are shown.



Device configuration



If your device has custom user interface (defined in the Device plugin script), it will be shown in the configure devices panel.



Pre-sampling warnings

When a user selects *Start sampling*, your plugin can temporarily halt the start sampling process to display a message.

In the example below, a plugin is telling the user that their device may be damaged if not used correctly. The user has the option of continuing (*Continue Sampling*), or canceling sampling.



Your hardware device.

Your hardware device provides one or more data streams to the host computer. The data streams may arrive in Lightning via USB, Bluetooth, serial, or other connection type.

Your plugin.

Your plugin is responsible for presenting the device streams to Lightning. It can allow users to configure your device and device streams through a graphical user interface that is defined in the plugin itself.

Lightning shows your device to the user in the devices list.

Lightning maps the data streams to signals and exposes these signals to users. Users are shown a graphical representation of your hardware device in the devices list.

The signals are displayed in Lightning channels.

